




Life Sciences

USD 2253a



# Palltronic® Integrity Test Instruments Spare Parts

# Palltronic Integrity Test Instruments – Spares List



**1½ in. Sanitary Cap with ¼ in. FNPT Threads\*** - Part #: FIT0023  
Used with: Palltronic Flowstar unit and Palltronic AquaWIT unit.



**6 mm ID x 8 mm OD Polyurethane Hose (Blue)** - Part #: S26.5.07  
Used with: Palltronic Flowstar unit and Palltronic AquaWIT unit.

*Note: Sold by the meter.*



**6 mm ID x 8 mm OD Polyurethane Hose (Red)** - Part #: S26.5.10  
Used with: Palltronic Flowstar unit and Palltronic AquaWIT unit.

*Note: Sold by the meter.*



**6 mm ID x 8 mm OD PTFE Hose “Steamable” (White)** -  
Part #: 64386  
Used with: Palltronic AquaWIT unit.

*Note: 9 meters standard.*



**AC Power Cord (110V U.S.)\*** - Part #: WIR0058  
Used with: Palltronic Flowstar unit and Palltronic AquaWIT unit.




**Data Transfer Cable, DB9** - Part #: S26.5.08  
Used with: Palltronic Flowstar unit and Palltronic AquaWIT unit.



**External Vent Valve Assembly** - Part #: S21.2.42  
Used with: Palltronic Flowstar unit and Palltronic AquaWIT unit.

*Note: Standard cord length is 1.5 meters. Longer cord lengths up to 20 meters are available. Please contact Pall for details*

# Palltronic Integrity Test Instruments – Spares List



**Stäubli Hose Fitting - Male (for 6 mm x 8 mm hose) - Part #: FIT0008**

Used with: Palltronic Flowstar unit and Palltronic AquaWIT unit.



**1 1/2 in. Sanitary Cap with Stäubli Male Adaptor - Part #: ACS0021BA**

Used with: Palltronic Flowstar unit and Palltronic AquaWIT unit.



**Stäubli Hose Fitting - Female (for 6 mm x 8 mm hose) - Part #: FIT0009**

Used with: Palltronic Flowstar unit and Palltronic AquaWIT unit.



**Paper Roll Spindle - Part #: S26.5.03**

Used with: Palltronic Flowstar unit and Palltronic AquaWIT unit.



**Stäubli Male Adaptor with 1/4 in. MNPT Threads\* - Part #: FIT0016**

Used with: Palltronic Flowstar unit and Palltronic AquaWIT unit.



**Suitcase with Interior Foam Padding**

Part #: S26.5.09a  
Used with: Palltronic Flowstar unit.



**Stäubli Female Adaptor with 1/4 in. MNPT Threads\* - Part #: FIT0017**

Used with: Palltronic Flowstar unit and Palltronic AquaWIT unit.



**External Vent Valve Assembly for Palltronic Flowstar MUX -**

*(Three types available)* – Used with: Palltronic Flowstar MUX unit.

Part # S29.5.01 – 3 meter length  
*(Supplied with MUX as standard)*

Part # S29.5.02 – 10 meter length

Part # S29.5.04 – 20 meter length



**Touch Screen Protector -**

Part #: S26.5.14

Used with: Palltronic Flowstar unit and Palltronic AquaWIT unit.



**Thermal Printer Paper Roll -**

*(Two types available)* – Used with: Palltronic Flowstar unit and Palltronic AquaWIT unit.

Part # S26.5.02 Standard paper – 24 meter length

Part # S26.5.02a Anti-static paper – 18 meter length



**Palltronic Flow Check Device -**

*(Three types available)* – Used with: Palltronic Flowstar unit and Palltronic AquaWIT unit.

Part # FC01  
(Flow Range: 8 mL/min ± 1)

Part # FC01M  
(Flow Range: 90 mL/min ± 10)

Part # FC01H  
(Flow Range: 900 mL/min ± 50)

\* These products are only available in the Americas



## Life Sciences

### New York - USA

+1 516 484 5400 phone

+1 516 801 9548 fax

biotech@pall.com e-mail

### Portsmouth - Europe

+44 (0)23 9230 3303 phone

+44 (0)23 9230 2506 fax

BioPharmUK@europe.pall.com e-mail

*Filtration. Separation. Solution.<sup>sm</sup>*

---


**Visit us on the web at [www.pall.com/biopharmaceutical](http://www.pall.com/biopharmaceutical)**

---

**Pall Corporation has offices and plants throughout the world in locations including:** Argentina, Australia, Austria, Belgium, Brazil, Canada, China, France, Germany, India, Indonesia, Ireland, Italy, Japan, Korea, Malaysia, Mexico, the Netherlands, New Zealand, Norway, Poland, Puerto Rico, Russia, Singapore, South Africa, Spain, Sweden, Switzerland, Taiwan, Thailand, United Kingdom, the United States and Venezuela. Distributors are located in all major industrial areas of the world.

Because of developments in technology these data or procedures may be subject to change. Consequently we advise users to review their continuing validity annually.

© 2008, Pall Europe Limited. Filtration. Separation. Solution. is a service mark of Pall Corporation.

 , Pall and Palltronic are trade marks of Pall Corporation. Stäubli is a trademark of Stäubli AG.

Part Numbers quoted above are protected by the Copyright of Pall Europe Limited.

® indicates a trademark registered in the USA, and ™ indicates a common law trademark.

PELWFR/SM/08.1830/01.2008